

Evaluating RF and HAMA Interference Recommended Testing Protocol

ASSURANCE™ RF and ASSURANCE™ HAMA from Sun Diagnostics, LLC are designed to simplify the evaluation of assays for interference from Rheumatoid Factor (RF) and Human Anti-Mouse Antibodies (HAMA). This recommended testing protocol is intended to guide the user through their RF and/or HAMA interference testing and assist them in interpreting their data.

1. For each interferent being evaluated, prepare three samples for testing:
 - a. Prepare 1:1 dilution of Interferent (RF or HAMA) and diluent.
 - b. Prepare a 1:1 dilution of sample and diluent. This is your CONTROL sample.
 - c. Prepare a 1:1 dilution of sample and interferent (RF or HAMA). This is your TEST sample.
2. Assay:
 - a. Assay each of the prepared samples using the method(s) or assay(s) being evaluated.
3. Evaluate:
 - a. Enter your appropriate data into the accompanying spreadsheet. We recommended using $\frac{1}{2}$ of the Total Allowable Error (TEa) as the Dmax (allowable difference allowed)
 - b. Review the output of the calculations.
 - c. The expected value is the mean of the CONTROL Sample (1b).
 - d. The Test result needs to be adjusted for any endogenous analyte present in the RF or HAMA material. The adjusted test result is calculated as 1c-1a.

For Technical Assistance or Support –

Email: support@sundiagnostics.us

Call: 1-877-SUN-DIAG (1-877-786-3424).