



LipoSep IP™ is a general laboratory reagent for isolation of HDL particles in human serum. LipoSep IP™ is for Laboratory Use only.

KIT COMPONENTS

Each Kit consists of 1 x 25 mL of precipitating reagent.

REAGENTS:

Immunoprecipitation reagent specific to human Apolipoprotein B.

NONREACTIVE INGREDIENTS

0.1% Sodium Azide

PRECAUTIONS AND WARNINGS

WARNING: Sodium Azide

Sodium Azide may react with lead and copper plumbing to form potentially explosive metal azides. Upon disposal, flush with a large volume of water to prevent azide build-up.

Dispose of all waste material in accordance with requirements of your local waste management authorities.

STORAGE AND STABILITY

This product will be stable until the expiration date when stored at 2 to 8°C. After use, promptly return to proper storage condition.

MATERIALS PROVIDED

LipoSep IP Immunoprecipitation reagent

MATERIALS REQUIRED BUT NOT PROVIDED

Pipettes, tips, and other lab consumables
Microcentrifuge, tubes

RECOMMENDED PROCEDURE:

1. Allow reagent to come to room temperature.
2. Pipette 200 µL of serum sample into a microcentrifuge tube.
3. Add 200 µL of LipoSep IP™ to the tube and vortex.
4. Incubate for 10 minutes at room temperature.
5. Centrifuge for 10 mins at 12,000 RPM in a bench-top microcentrifuge.
6. After centrifugation, carefully remove the supernatant from the tube. Do not disturb the pellet at the bottom.
7. The supernatant contains the HDL particles from the sample.
8. Assay as desired. Any results should be multiplied by 2 to obtain value of original sample.

NOTE: A larger sample volume may be used as long as an equal volume of LipoSep is also added.

LIMITATIONS

- This product should not be used past the expiration date listed on the kit box and storage container.
- If there is evidence of microbial contamination discard the bottle.
- This product is for laboratory use only. Not for drug, household or other use.

For customer or technical support:

email support@sundiagnosics.us

phone: 877-SUN-DIAG (877-786-3424)



Sun Diagnostics, LLC.
60 Pineland Drive
Brunswick Hall, Suite 322
New Gloucester, ME 04260
Toll free: 1-877-786-3424
www.sundiagnosics.us

Symbols Used:			
	Lot Number		Consult Package Insert
	Expiration		Manufacturer
	Storage Temp		Reorder Number
LipoSep™ IP is a general purpose reagent intended for laboratory use. CE marking not required			