

ASSURANCE™ Interference Test Kit is for Laboratory Use as part of Interference testing procedures. Interferents available include: CREATININE, GLUCOSE, GLYCEROL LACTATE, PYRUVATE, BLOOD UREA NITROGEN (BUN), URIC ACID and 0.1N NaOH as the base for URIC ACID. Components may be purchased individually, in customer-defined combinations, or as INT-06 (one vial of each component).

KIT COMPONENTS AND TYPICAL VALUES

Interferent	Typical Concentration, mg/dL	Reorder Number
CREATININE	300	INT-06C
GLUCOSE	20,000	INT-06GLU
GLYCEROL	340	INT-06GLY
LACTATE	1,800	INT-06LAC
PYRUVATE	100	INT-06PYR
BLOOD UREA NITROGEN	2,400	INT-06BUN
URIC ACID	470	INT-06UA
NaOH <i>(for UA control pool)</i>	0.1N	
ASSURANCE™ Interference Test Kit – 1 Vial of each Component		INT-06

REAGENTS

This product is prepared from purified chemicals. Human source materials are not used.

NONREACTIVE INGREDIENTS

None

PRECAUTIONS AND WARNINGS

Dispose of all waste material in accordance with requirements of your local waste management authorities.

STORAGE AND STABILITY









This product will be stable until the expiration date when stored at -10 to -20°C. After use, refreeze promptly. A maximum of 4 Freeze/Thaw cycles is recommended.

LIMITATIONS

- Typical concentrations are intended as a guide only. **Laboratory should determine actual values before proceeding with interference experiment(s).**
- INT-06 should not be used past the expiration date listed on the kit box. Refer to Interferent vials for expiration of individually purchased components or custom kit configurations.
- If there is evidence of microbial contamination discard the vial.
- This product is not intended for use as a control, standard or calibrator.
- Information gathered from experiments performed using ASSURANCE™ Interference Test Kit should not be used to adjust patient

results. Confirm results with an alternate method when necessary.

Sample Interference Testing procedure (based on CLSI EP-7A guidelines), dilution calculator and results reporting spreadsheet are available at www.sundiagnosics.us or email support@sundiagnosics.us.

Symbols Used			
	Biohazard		Lot Number
	Protect from light		Consult Package Insert
	Manufacturer		Reorder Number
	Expiration	ASSURANCE™ Interference Test Kit is a general purpose reagent intended for laboratory use. CE marking not required	
	Storage Temp		