



Recent QC Materials:

| Analyte             | Matrix   | Analyte Source                       |
|---------------------|--|--------------------------------------|
| HDL <sub>3</sub> -C | Human serum-based with preservatives and stabilizers | Purified human lipoproteins          |
| sdLDL-C             | Human serum-based with preservatives and stabilizers | Purified human lipoproteins          |
| MPO                 | Human serum-based with preservatives and stabilizers | Purified MPO from human neutrophils  |
| Leptin              | Human serum-based with preservatives and stabilizers | Recombinant human                    |
| Lp-PLA <sub>2</sub> | Human serum-based with preservatives and stabilizers | Human serum                          |
| GAL-3               | Human serum-based with preservatives and stabilizers | Recombinant human                    |
| Lp(a)               | Human serum-based with preservatives and stabilizers | Purified human Lp(a)                 |
| Triglycerides       | Human serum-based with preservatives and stabilizers | Human triglyceride-rich lipoproteins |
| Glycerol            | Human serum-based with preservatives and stabilizers | Glycerol                             |

Contact [support@sundiagnosics.us](mailto:support@sundiagnosics.us) for more information about these, or any other custom QC needs.