



## Technical Bulletin #19-005

<b>Date</b>	May 10, 2019
<b>Part Numbers</b>	INT-04 and INT-06
<b>Product Name</b>	ASSURANCE™ Drug Interference Test Kit and components ASSURANCE™ Endogenous Interference Test Kit and components
<b>Subject</b>	Updated Concentrations

Dear Valued Customer,

In order to comply with new recommendations from the Clinical Laboratory Standards Institute (CLSI) we have changed many of the concentrations of interference materials in the INT-04 and INT-06 ASSURANCE™ Interference Test Kits.

The new concentrations allow a convenient 20-fold dilution of interference material to reach the recommended CLSI test concentration.

Please consult the package inserts for concentration information and contact Sun Diagnostics, LLC with any questions.

Regards,

Technical Support  
Sun Diagnostics, LLC  
Email: [support@sundiagnosics.us](mailto:support@sundiagnosics.us)  
Phone: 1-877-786-3424

TB19-005 INT-04.06 UPDATED CONCENTRATIONS