

**ASSURANCE™ Triglyceride-Rich Lipoproteins** simplify interference testing and more accurately reflect lipemia. Our pooled human TRLs contain a range of chylomicron, VLDL, IDL, and LDL particle sizes just as you would see in human sera.

Triglycerides	>15,000	mg/dL
CM+VLDL, Total	8493	nmol/L
CM+Large VLDL	822	nmol/L
Medium VLDL	4775	nmol/L
Small VLDL	2895	nmol/L
LDL, Total	25,500	nmol/L
IDL	5250	nmol/L
Large LDL	1083	nmol/L
Small LDL	19,167	nmol/L
Average VLDL Size	55.7	nm
Average LDL Size	19.9	nm

